Instrument: Pegasus® GC-HRT<sup>+</sup>

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## 5. References

<sup>1)</sup>Anstett A., Fanny C., Alonso D.E., and Waddell Smith R., "Characterization of 2C-phenethylamines using high-resolution mass spectrometry and Kendrick mass defect filters", *Forensic Chemistry*, **2018**, 7, 47-55.

<sup>2)</sup>Dubois L.M., Pierre-Hugues S., Heudt L., Focant J.F., and Perrault K.A., "Characterizing decomposition odor from soil and adipocere samples at death scene using HS-SPME-GCxGC-HRTOFMS", *Forensic Chemistry*, **2018**, 8, 11-20.

<sup>3)</sup>Veenas C., Ripszam M., Glas B., Liljelind I., Claeson A.S., and Haglund P., "Differences in chemical composition of indoor air to rooms associated/not associated with building related symptoms", *Science of the Total Environment*, **2020**; <https://doi.org/10.1016/j.scitotenv.2020.137444>.

<sup>4)</sup>Mazur D.M., Detenchuk E.A., Sosnova A.A., Artaev V.B., and Lebedev A.T., "GC-HRMS with complementary ionization techniques for target and non-target screening for chemical exposure: Expanding the insights of the air pollution markers in Moscow snow", *Science of the Total Environment*, **2021**, 761, 144506-144517.

<sup>5)</sup>Lebedev A.T., Mazur D.M., Artaev V.B., and Tikhonov G.Y., "Better screening of non-target pollutants in complex samples using advanced chromatographic and mass spectrometric techniques", *Environmental Chemistry Letters*, **2020**, 18, 1753-1760.

<sup>6)</sup>Daygon V.D., Calingacion M., Forster L.C., De Voss J.J., Schwartz B.D., Ovenden B., Alonso D.E., McCouch S.R., Garson M.J., and Fitzgerald M.A., "Metabolomics and genomics combine to unravel the pathway for the presence of fragrance in rice", *Scientific Reports*, **2017**, 8, 8767; DOI:10.1038/s41598-017-07693-9

<sup>7)</sup>Di Giovanni N., Meuwis M.A., Louis E., and Focant J.F., "Untargeted serum metabolic profiling by comprehensive two-dimensional gas chromatography-high-resolution time-of-flight mass spectrometry", *J. Proteome Res.*, **2020**, 19(3), 1013-1028.

<sup>8)</sup>Ncube S., Dube S., and Nindi M.M., "Determination of volatile compounds during deterioration of African opaque beer using a stir bar sorptive extraction technique and gas chromatography-high resolution mass spectrometry", *Current Research in Food Science*, **2020**, 3, 256-267.

<sup>9)</sup>Vyviurska O., and Spanik I., "Assessment of Tokaj varietal wines with comprehensive two-dimensional gas chromatography coupled to high resolution mass spectrometry", *Microchemical Journal*, **2020**, 152, 104385; <https://doi.org/10.1016/j.microc.2019.104385>.

<sup>10)</sup>Alisi I., Zoccali M., Tranchida P.Q., and Mondello L., "Analysis of organic compounds in coal tar by using comprehensive two-dimensional gas chromatography-high resolution time-of-flight mass spectrometry", *Separations*, **2020**, 7, 26; <https://doi.org/10.3390/separations7020026>.

<sup>11)</sup>Gross J.H., *Mass Spectrometry – A Textbook*, Springer, **2004**, Chp. 7.

**LECO Corporation** | 3000 Lakeview Avenue | St. Joseph, MI 49085 | Phone: 800-292-6141 | 269-985-5496  
info@leco.com • www.leco.com | ISO-9001:2015 Q-994 | LECO is a registered trademark of LECO Corporation.  
Pegasus, ChromaTOF are registered trademarks of LECO Corporation.