

REAGENT TUBES

- 771-030 Glass Tube, requires two (2) 772-520 O-rings for sealing
- 771-373 Glass Tube, requires two (2) 772-520 O-rings for sealing
- 618-603 Quartz Catalyst Tube, requires two (2) 772-520 O-rings for sealing, and a 772-738 O-ring for the removable block
- 772-739 Glass Catalyst Tube, requires two (2) 772-520 O-rings for sealing, and a 772-738 O-ring seals the removable block

OPERATION ACCESSORIES

- 619-591-877 Balston Filter
- 761-929 Crucible Tongs
- 789-149 EMI Gasket
- 775-307 Filter Removal Tool (Thumbscrew)
- 775-306 Secondary Filter, 10 micron
- 503-032 Glass Accelerator Scoop
- 773-579 Metal Accelerator Scoop
- 764-078 One gram Weight
- 760-138 Sample Tweezers
- 619-825 Trap Tube for CSLS600 (only)
- 501-241 Vacuum Grease

*NOTE: Part numbers and standards' values may change.
Consult LECO for latest information.*