


RH/ROH/TCH600-Series Consumables & Spare Parts Reference Card

OBSOLETE—No further updates will be available.



See reverse side for contact information

FUSION SUPPLIES

501-598 Etchable Nibbled Nickel 100 g bottle
501-073 Graphite Powder 100 g bottle
618-376* High Temperature Electrode Insert (required for
782-720 Crucibles)
502-040 Large Tin Capsule 100/bottle
774-693 Lower Electrode Assembly, requires 775-312 O-ring,
includes: 617-866 Electrode (617-866 for EN models),
611-351-181 Electrode Insert, 774-690 Knurled Ring (Collar)
502-822 Nickel Capsule 100/bottle
763-029 1.5 g Nickel Baskets 100/bottle
502-345 1.5 g Nickel Baskets Ultra High Purity 100/bottle
502-344 1.0 g Nickel Baskets Ultra High Purity 100/bottle
501-059 Small Tin Capsule 100/bottle

**Not part of 774-693 Assembly*

CATALYSTS AND REAGENTS

501-171-HAZ Anhydron 454 g (1 lb.) bottle
502-705 Copper Sticks 100 g ampule
502-878 Copper Sticks 100 g bottle
502-656 Copper Turnings 60 g bottle
502-190 Copper Oxide Wire 100 g bottle
501-081 Glass Wool 454 g (1 lb.) box
501-073 Graphite Powder 100 g bottle
502-174-HAZ Lecosorb 500 g bottle, 20-30 mesh
502-176-HAZ Lecosorb 500 g bottle, 8-20 mesh

OPERATION ACCESSORIES

766-053 Crucible Tweezers
766-646 Air Filter Coating
775-306 Secondary Filter, 10 micron
775-307 Screen Filter Removal Tool (Thumbscrew)
625-711-599 Coolant Additive
764-078 One gram Weight
783-785-110-HAZ Replacement OMI Tube
760-138 Sample Tweezers (Straight)
501-241 Vacuum Grease
767-473 Utility Funnel
620-760 Vacuum Filter (Ash Vac)
618-809 Medium Reagent Tube (Glass)

*NOTE: Part numbers may change.
Consult LECO for latest information.*

*Please Note: Quartz Tubes are typically utilized for higher
temperature applications (i.e. catalyst tubes). Glass reagent/filter tubes
are designed for lower temperature applications. Glass tubes should
not be used in place of quartz tubes.*