



**FUSION SUPPLIES**

501-598 Etchable Nibbled Nickel 100 g bottle  
501-073 Graphite Powder 100 g bottle  
502-040 Large Tin Capsule 100/bottle  
611-351-180 Lower Electrode Tip (for 764-330 Crucibles with Automation)  
611-351-183 Lower Electrode Tip (for 764-330 Crucibles with Automation)  
619-896 Lower Electrode Tip (for 619-895 Crucibles)  
502-822 Nickel Capsule 100/bottle  
502-345 1.5 g Nickel Baskets Ultra High Purity 100/bottle  
502-344 1.0 g Nickel Baskets Ultra High Purity 100/bottle  
501-059 Small Tin Capsule 100/bottle  
761-744 Molecular Sieve Assembly

**CATALYSTS AND REAGENTS**

501-171-HAZ Anhydrous 454 g (1 lb.) bottle  
502-705 Copper Sticks 100 g ampule  
502-878 Copper Sticks 100 g bottle  
502-995 Copper Turnings 60 g bottle  
501-081 Glass Wool 454 g (1 lb.) box  
501-073 Graphite Powder 100 g bottle  
502-174-HAZ Lecosorb 500 g bottle, 20-30 mesh  
502-176-HAZ Lecosorb 500 g bottle, 8-20 mesh  
761-747-HAZ Schutze Reagent

**OPERATION ACCESSORIES**

766-053 Crucible Tweezers  
766-646 Air Filter Coating  
775-306 Secondary Filter, 10 micron  
775-307 Screen Filter Removal Tool (Thumbscrew)  
625-711-599 Lubricool Cooling Agent  
764-078 One gram Weight  
783-785-110-HAZ Replacement Tube for 783-786 Kit  
760-138 Sample Tweezers (Straight)  
501-241 Vacuum Grease  
767-473 Utility Funnel  
620-760 Vacuum Filter (Ash Vac)  
618-809 Medium Reagent Tube (Glass)

*NOTE: Part numbers and standards' values may change. Consult LECO for latest information.*

*Please Note: Quartz Tubes are typically utilized for higher temperature applications (i.e. catalyst tubes). Glass reagent/filter tubes are designed for lower temperature applications. Glass tubes should not be used in place of quartz tubes.*